

**ViewPoint**V6™  
TECHNOLOGY WALL STATION

**Revolutionizing the way caregivers interact  
with technology ... and their patients.**



Humanscale™  
*Healthcare*

An important addition to Humanscale Healthcare's ViewPoint line, the V6 technology wall station combines exceptional space-savings with leading ergonomic design for comfortable caregiver computer use while sitting or standing. The V6 supports a monitor and keyboard—along with hardware, shelves and other accessories—within an ergonomic adjustment range for easy positioning to provide a healthy computing experience and encourage engagement with patients. Its sealed cable management system also promotes organization and mitigates infection control concerns.



## Features

- Flexible monitor arm options facilitate and improve patient-caregiver interactions
- Wide range of dynamic adjustability accommodates most users in sitting and standing positions
- Ergonomic design encourages comfortable and convenient computer use
- Numerous monitor and keyboard arm lengths for ultimate versatility
- Sophisticated and customizable cable management system offers easier cleaning, reduces infection control concerns and prevents the ingress of liquids
- Accessory mount supports CPU holder or rail, which can accommodate peripheral shelves, and mounts for gloves, wipes and other common items
- Streamlined aesthetic features anodized aluminum and a sophisticated color palette, complementing all environments
- Sustainable composition allows unit to contribute to the achievement of valuable LEED credits

## Specifications

- More than 11 inches of near effortless monitor and keyboard height adjustability
- Customizable track is available in five standard lengths: 27 inches, 37 inches, 47 inches, 57 inches or 67 inches
- Dynamic monitor arm supports single or dual monitors up to 40 pounds total
- 180° stop mechanism on arms prevents wall damage
- Monitor extends up to 26.25 inches from wall
- Keyboard (including platform) extends up to 36.75 inches from wall
- Folds to a maximum depth of 5 inches from wall
- Innovative 15° negative tilt keyboard support encourages ergonomic use of standard or custom keyboards
- 5-year warranty